



# BRIT Systems DoctorWorks iPad Viewer



*Designed specifically for the iPad*

*View:*

*Patient Lists*

*Histories*

*Reports*

*Studies*

*Urgent Findings*

*ER Calls*

DoctorWorks iPad Viewer connects with a Roentgen Works Server to provide healthcare workers portable access to radiology studies and reports. Users can also participate in Urgent Findings and ER Discordances according to their role: radiologists can enter urgent findings and respond to ER preliminary diagnoses. Other healthcare professionals can view and acknowledge Urgent Findings; enter ER preliminary diagnosis and acknowledge the radiologist's discordance.




DoctorWorks for the iPad was designed specifically for this device, optimizing the screen space and utilizing the iPad gestures for navigation and image viewing. The application runs locally on the iPad device, utilizing the Roentgen Works servers for authentication, HIPAA logging and access to other servers. Roentgen Works provides the DICOM hosting and integration tools and HIE Patient Discovery and MPI tools to give the iPad access to all studies within an enterprise and cross enterprises.

After authenticating, users view their default worklist as established in Roentgen Works. There is also a **My Patients** worklist for quick access to all patients listed as having recent studies ordered by that physician. Studies with unacknowledged Urgent Findings and ER Discordances display icons on the worklist. Once a study has been selected, the Urgent Finding / ER Discordance details appear on the screen.


**Custom Solutions for Medical Imaging**

1909 Hi Line Dr., Dallas, TX 75207 (800) 230-7227 Phone: (214) 630-0636 www.brit.com

# DoctorWorks iPad Viewer



Open studies via clicking on the image icons.



## Image Viewing Tools Include:

- Auto-synchronization of series
- Synchronized roaming of series
- Image magnification via spreading and closing fingers on screen
- Roaming via single finger drag
- W/L via selection of presets and single finger drag

## ER Discordance and Urgent Finding Tools:

- Enter Urgent Findings
- Enter ER Calls and Discordances
- Acknowledge Urgent Findings and Discordances
- Review Histories of events



## Extensive Patient/Study Search Options:

- Name, ID, Accession #
- Description, CPT
- Modality Type
- Body Region
- Study Status
- Study Date Range
- Updated Date Range
- Patient Age Range

System Requirements: Runs on any iPad. Requires Roentgen Works Server

DoctorWorks and Roentgen Works are Trademarks of BRIT Systems

**Custom Solutions for Medical Imaging**

1909 Hi Line Dr., Dallas, TX 75207 (800) 230-7227 Phone: (214) 630-0636 www.brit.com